
XDCOLLECTION

**Revix RCS recycled
plastic massage gun**

Content

English	_____	p. 3
Nederlands	_____	p. 6
Deutsch	_____	p. 9
Français	_____	p. 12
Español	_____	p. 15
Svenska	_____	p. 18
Italiano	_____	p. 21
Polski	_____	p. 24
EU Declaration of conformity	_____	p. 27
Explanation WEEE logo	_____	p. 28

ENGLISH

DEVICE LAYOUT:

1. Power button
2. Handle
3. LED indicator
4. Type-C Input
5. Massage head

SPECIFICATIONS:

Type-C input	5V/1A
Battery capacity	1200mAh / 3.7V
Charging time	About 2.5 hours
Massage type	Vibrate
Massage amplitude	6mm
Max. Movements per minute	8000
Min. Movements per minute	4000
Max. Noise level	60dB
Min. Noise level	50dB
Run time	3.5 hours

*Times are approximate times only, depending on battery condition.

Please read the following instructions carefully and follow them to ensure that the product serves you safely and to your full satisfaction.

IMPORTANT:

1. Before first use first use, fully charge your massage gun until the red LED indicator light changes to blue.
2. Only use included accessories or certified and tested accessories to power this device. Using third party accessories can influence the item safety.

INSTRUCTIONS:**Charging:**

Connecting to Power: Insert the connector of the included USB-C power cable into the power input port on the back of the product. Insert the other connector into the USB-C port on a computer or into a compatible 5V/1A USB-C power adapter (not included) and plug it into a wall outlet. While charging the massage gun, the red LED will light up. Once the massage gun is fully charged, the LED light will turn blue.

Power button:

- Long press power button (⏻) to turn the massage gun ON. Long press again to turn OFF the massage gun.
- Powering on the device automatically sets the intensity level to 1. Short Press (⏻) to change the intensity level. There are 2 intensity levels.

SAFETY PRECAUTIONS:

When using this product, basic precautions should always be followed:
Read all the instructions before use.

WARNING TO REDUCE THE RISK OF BURNS, FIRE, ELECTRIC SHOCK OR INJURY TO PERSONS:

- Always unplug this device from any electrical outlets before using and before cleaning.
- An electrical device should never be left unattended when plugged in. Unplug from the outlet when not charging.
- Do not operate under a blanket or pillow. Excessive heating can occur and cause the fire, electric shock, or injury to persons.


- Use this device only for its intended use as described in this manual.
- Never operate this device if it is damaged, not working properly, if it has been dropped or has come in contact with water.
- Do not operate where aerosol (spray) products are being used or where oxygen is being administered.

- When charging, connect this device to a properly grounded outlet.
- Intended use of this personal massage gun is for gentle massage only.
- Avoid applying excessive force to the massage gun. This can result in overload and can cause damage to the massage gun and/or to the user.
- Do not use this massage gun in close proximity to loose clothing, jewelry, hair or anything that can get caught or pinched in the massage gun.
- Keep away from children, this is NOT a toy.
- Do not use this massage gun while it is being charged.
- Avoid any areas of skin/body that may be injured.
- Always use caution when operating this massage gun to avoid injury.

DANGER

- Do not reach for an electrical device that has fallen into water. If this occurs while the device is plugged in, unplug immediately.
- Do not use while bathing or in a shower.
- Do not place or store where it can fall or be pulled into a tube or sink. Do not place in or drop into water or other liquid.

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EU Declaration of Conformity

ENGLISH

Xindao B.V. hereby declares that this product is in compliance with the essential requirements and other relevant provisions of Directive 2014/30/EU.
The complete Declaration of Conformity can be found at www.xindao.com > search on item number.

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Explanation WEEE logo:

ENGLISH

The crossed-out wheeled bin symbol indicates that the item should be disposed of separately from household waste. The item should be handed in for recycling in accordance with local environmental regulations for waste disposal. By separating a marked item from household waste, you will help reduce the volume of waste sent to incinerators or land-fill and minimize any potential negative impact on human health and the environment.

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XD Connects B.V.